

# SKINCARE ROUTINE IMASK EMS HANDSET

## -ALL SKINTYPES-



### DAILY CARE

Gentle Face Wash



Blue Photon  
Daily  
Toning & Energise  
Collagen Synthesis  
Mode with skin  
toner cotton pad



Apply other skincare  
serum/ moisturizer as norm

### DAILY CARE




Blue Photon  
Daily or As Required

Makeup Remove & Gentle  
Exfoliate  
Clean Mode with Cleansing  
Water or Makeup Remover  
cotton pad



Recommended  
for all skin  
types Daily

Gentle Face Wash  
Extra fine brush knot clears  
dirt & pollution collected in  
pores more effectively &  
efficiently



### WEEKLY COMPLETE REJUVENATION TREATMENTS



Blue or Red Photon  
1-2 times per week

Skin Care Preparation Pore  
Relaxation & Opening  
Moisturizing Mode  
Apply Facial Toner and  
massage in circular motion




Highly Recommended  
for Oily &  
Combination Skin  
1-2 times per week

iMask Microscope Pore  
Cleaner  
Add target pore treatment to  
clear pollution and unblock  
pores. Monitor and camera  
record skin condition on your  
mobile device



Red Photon  
1-2 times per week

Deep Rejuvenation  
EMS Lifting Mode  
Apply Facial Mask or Facial  
Massage Oil to face & neck.  
Massage in circular motion



Red Photon  
Daily

Skin Care Import  
Moisturizing Mode  
Apply Facial  
serum/moisturizer and  
massage in circular motion



Blue or Red Photon  
Daily

Balancing & Pore Shrink  
Cool Mode  
Massage by tapping over  
face to lock-in nutrient & skin  
moisture

